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
▲instruments

☎ 802 431 408

www.controltecnica.com

Manuale d'uso
Instruction Handbook

Congelatori
Freezers

PLATINUM

110 SH - 110 H
370 SH - 370 H
550 SH - 550 H

340 SV - 340 SV Plus
340 V - 340 V Plus

500 SV - 500 SV Plus
500 V - 500 SV Plus

Cod. 514903 Rel. 171207

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▲instruments

☎ 802 431 408

www.controltecnica.com

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IMPORTANT WARNING

The high level of technology of the products made by the manufacturer occasionally requires components to be adapted from parts normally found on the market in their standard commercial version. Therefore, if it should be ascertained that any parts have been replaced with other components that are not original, i.e. that have not been supplied and guaranteed directly by the MANUFACTURER, or that any unauthorised technical interventions have been carried out, the guarantee, if any is in force, shall be considered no longer valid and the MANUFACTURER himself shall no longer be held responsible towards the purchaser and towards third parties.

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The machine complies with the regulations in the following Italian and European Directives:

- D.P.R. 459/96 regulation for implementing Directives 89/392/EEC, 91/368/EEC, 93/44/EEC and 93/68/EEC concerning the approximation of the laws of the Member States relating to machines.
- Machinery Directive n° 98/37/EEC (22-06-1998) and subsequent modifications.
- Low-Voltage Directive 2006/95/EC.
- Electromagnetic Compatibility Directive n° 2004/108/EC.

In particular it complies with the following Harmonized Standards:

UNI EN ISO 12100-1:2005: Safety of machinery - Basic concepts, general principles for design - Part 1: Basic terminology, methodology.

UNI EN ISO 12100-2:2005: Safety of machinery - Basic concepts, general principles for design - Part 2: Technical principles.

UNI EN 294:1993: Safety of machinery. Safety distances to prevent danger zones being reached by the upper limbs.

CEI EN 60204-1: Safety of machinery. Electrical equipment of machines. Part 1: General requirements. ed. 1998

CEI EN 61010-1: Safety requirements for electrical equipment for measurement, control and laboratory use - Part 1: General requirements. ed. 2001

CEI EN 61326-1: Electrical equipment for measurement, control and laboratory use - EMC requirements - Part. 1 - General requirements. ed. 1998 and subsequent modifications A1 (1999) A2 (2002) A3 (2004).

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The materials supplied are certified:

- to conform in every way with the relative drawings, specifications and order;
- to have been checked and tested in conformity with the conditions and requirements of the order, according to international safety standards;
- pursuant to the provisions of Directive 93/42/EEC concerning medical devices, even though the equipment is not to be considered a medical device in the strict sense, as it is not intended for use in direct contact with the patient.

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SUMMARY

1	SUMMARY OF RATING PLKATE DATA.....	2
2	WARNINGS.....	3
2.1	General warnings.....	3
2.2	Warnings for transport and handling.....	4
2.3	Warnings for installation.....	4
2.4	Warnings for personnel in charge of the machine.....	4
2.5	Warnings for maintenance.....	5
2.6	Explanation of symbols.....	5
3	TECHNICAL SPECIFICATIONS.....	7
3.1	Technical data.....	7
3.2	Fuses.....	8
3.3	Environmental conditions.....	8
3.4	Technical lay-out diagrams.....	9
4	HANDLING AND REMOVAL OF PACKAGING.....	10
4.1	Personnel requisites.....	10
4.2	Machine conditions.....	10
4.3	Equipment needed for handling.....	10
4.4	How to raise and move with a bridge crane.....	11
4.5	Removal of packaging.....	11
4.5.1	How to remove the outer straps.....	12
4.5.2	How to remove the cardboard casing.....	12
4.5.3	Removal of plastic materials.....	12
4.5.4	How to remove the wooden side blocks.....	12
4.5.5	How to remove the pallet.....	13
4.6	How to block the machine.....	13
5	DESCRIPTION OF THE SYSTEM.....	15
5.1	General view.....	15
5.2	Internal view.....	16
5.3	Description of external connections.....	17
5.4	Control system.....	18
5.5	"Bioguard" system.....	19
5.6	Cooling system.....	19
5.7	Functioning of the cooling system.....	19
5.8	Voltage booster (optional).....	20
5.9	Functioning of the auxiliary LN ₂ /CO ₂ cooling system (optional).....	20
5.10	Insulation.....	20
6	INSTALLATION.....	21
6.1	Positioning of the machine.....	21
6.2	How to connect the LN ₂ /CO ₂ cylinders (optional).....	23
6.2.1	How to replace the cylinder.....	23
6.3	Electrical wiring.....	24
6.3.1	How to connect to the electrical mains supply.....	24
6.3.2	Replacement of the supply cable.....	25
6.3.3	Remote alarm connection.....	25
6.3.4	Serial interface (optional).....	26
7	USE FORESEEN BY THE MANUFACTURER.....	27
7.1	Aim of the machines.....	27
7.2	Intended use of the machines.....	27
7.3	Operator.....	27
7.4	Use limits.....	27
7.5	Use of protective clothings.....	27
7.6	Dangerous areas and waste risks.....	27
8	START-UP.....	28
8.1	Description of the control panel.....	28
8.2	Start-up.....	30
8.3	Setting the operating temperature.....	30
8.4	Setting the alarm temperature set point.....	31
8.5	Activation/deactivation of the Bioguard® system and locking/unlocking the control panel (Only for models equipped with this device).....	32
8.6	Card memorisation (tag) - (Bioguard® models only).....	33
8.7	How to open/ close the door.....	34
8.7.1	Description of how to open/close the door with a key.....	34
8.7.2	Description of how to open/close the door with the Bioguard® system.....	34
9	USE.....	35
9.1	How to install the drawers (optional) - vertical models only.....	35
9.2	How to install the shelves- vertical models only.....	36
9.3	Product loading.....	36
9.4	Temperature chart recorder (optional).....	36
9.4.1	Use of the recorder.....	37
9.5	Data downloading via serial port (optional).....	39
9.6	Additional probe (optional).....	39
9.7	Port hole (optional).....	39
9.8	Water condenser connection (optional).....	39
9.9	Storage containers.....	40
10	SAFETY DEVICES- CHECK AND SET-UP.....	41
11	MAINTENANCE.....	42
11.1	How to defrost and clean the main body and doors.....	42
11.2	Cleaning of the air condenser.....	43
11.3	Cleaning of the water condenser.....	43
11.4	Cooling system maintenance.....	43
11.5	Maintenance of the electrical-electronic part.....	43
12	TROUBLESHOOTING.....	44
13	REMOVAL FROM INSTALLATION SITE.....	46
13.1	Disassembly.....	46
13.2	Scrapping.....	46
14	DIAGRAMS.....	47
14.1	Refrigerator drawing (-40 °C).....	47
14.2	Refrigerator drawing (-80 °C).....	48
14.3	Electric diagram (-40 °C).....	49
14.4	Electric diagram (-85 °C).....	51
14.5	Electric diagram (equipment at -85 °C with booster).....	53

  		2 Angelantoni INDUSTRIE MORE THAN TECHNOLOGY	
Model 3			
Serial No	4	Year	5 Weight kg 8
C	6	RH%	7 IP 9
POWER SUPPLY	V~	10	Hz 13
	P Kw	11	A 12 Class 14
REFRIGERANT PROPERTIES	n.1	15	Pmax kPa 17
	n.2	16	Pmax kPa 18
COUNTRY OF ORIGIN: ITALY			

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GB 1 SUMMARY OF RATING PLATE DATA

- Find the special rating plate on the machine in order to check the technical data.
- Check the machine model and supply voltage before carrying out any operation whatsoever.
- If you find any discrepancies, contact the manufacturer or your supplier immediately.

1	EC compliancy mark	10	Supply voltage
2	Manufacturer's address	11	Absorbed power
3	Model	12	Electrical absorption
4	Serial N°	13	Frequency
5	Year of manufacture	14	Isolation class
6	Temperature	15-16	Type of refrigerant
7	Relative humidity	17-18	Pressure of refrigerant
8	Weight		
9	Electrical protection grade		

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GB 2 WARNINGS

2.1 GENERAL WARNINGS

- Do not carry out any operations or manoeuvres unless you are absolutely certain of their effect; if in doubt, contact your nearest technical assistance service or the manufacturer himself.
- The manufacturer will not be held responsible for damage to the machine or to objects in the following cases:
 - improper use of the machine
 - use of unsuitable personnel
 - incorrect assembly and installation
 - defects in the plant systems
 - unauthorized modifications or operations to the machine
 - use of spare parts that are not original pieces
 - failure to comply with the norms given in this handbook
 - exceptional events
- This instruction handbook has been designed for the following personnel:
 - Personnel in charge of transport, handling and removal of packaging
 - Personnel in charge of the preparation of the plant systems and installation site
 - Installers
 - Personnel in charge of using the machine
 - Personnel in charge of maintenance
- The instruction handbook indicates the use foreseen by the manufacturer and cannot ever replace adequate experience of the operator. It can only be used as a reminder of the main operations to be carried out.
- The instruction handbook should be kept carefully and should also be within easy reach for reference. If necessary, photocopy the pages concerned directly with machine use.

The handbook should last at least the life-time of the machine itself.

- The instruction handbook gives technical information on how the machine is manufactured at the present time; the manufacturer reserves the right to carry out any modifications he deems necessary to the machines and to the instruction handbooks, without prior notice or replacement. If you lose or destroy the handbook, you may ask for a copy. Please give the rating plate data (see chap. 1) in your request.

This machine is not referred to in enclosure IV of the community norm on machines, and therefore the procedure in article 8, paragraph 2, letter A in the Community norms for machines 89/392/CEE and 91/368/EEC has been applied. Article 8, paragraph 2, letter A of the Community norms for machines 89/392/CEE and 91/368/EEC obliges the manufacturer to provide the file foreseen by the enclosure of the afore-mentioned laws, and to store it in the technical archives of the company Angelantoni Industrie, località Cimacolle - Massa Martana (Perugia).

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GB 2 WARNINGS

2.2 WARNINGS FOR TRANSPORT AND HANDLING


- This symbol, placed on each packaging, indicates the weight of each package.
- Always check that the tools and machines to handle and transport the machine are adequate.
- Always keep the machine in an upright position. If the machine should accidentally turn upside-down or on its side, do not switch it on. Put it in the correct position and contact the manufacturer.

2.3 WARNINGS FOR INSTALLATION

- Installation should always be carried out by specialized personnel.
- Carefully follow the instructions on how to prepare the plant systems before installing the machine.
- When the installation site is being prepared, bear in mind the space and work conditions of the personnel in charge of the machine so as to reduce to a minimum noise, fatigue, discomfort and anything else which may have a negative influence on the staff.
- When designing the installation site, remember to leave sufficient space for control, maintenance, cleaning, and removal of production waste material.
- Make sure that the work site is adequately lit so that personnel can work in optimum conditions.
- When designing the installation site, please refer to the norms in force and in particular:
 - set up all the firefighting and safety devices.

2.4 WARNINGS FOR PERSONNEL IN CHARGE OF THE MACHINE

- The machine may only be used by personnel who have read the rules described in this handbook.
- To open the chamber when its temperature is very different from the ambient temperature could cause problems: if the internal temperature is very low, condensate and frosting could be caused by ambient humidity. If the door is often opened under these conditions, the heat exchangers could be obstructed.

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GB 2 WARNINGS

2.5 WARNINGS FOR MAINTENANCE

- Always disconnect the machine from the electrical mains before carrying out any maintenance operation.
- Do not use solvents or alcohol to clean the varnished parts as these products could damage the surface.

2.6 EXPLANATION OF SYMBOLS

- The symbols shown below may be found on the machine or in this instruction handbook.
- Pay attention to their meaning before going any further.

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- 1 **Use gloves**
- 2 **GENERAL DANGER**
This symbol is always accompanied by an explanation of the danger
- 3 **WARNING**
Important note
- 4 **Operations that must be carried out by at least two people**

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
⁽¹⁾ Performance is given with the freezer: V - H series stabilised at -80 °C and SV - SH series stabilised at -40 °C (60068-3-5 IEC normative).

⁽²⁾ Measured in front of the source at a distance of 1 m, and a height of 1,3 m, in a non-reverberating environment.

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
3.3 ENVIRONMENTAL CONDITIONS
 In order for the appliance to operate correctly, it should be placed in a site with the following requisites:

- far from heat sources,
- far from direct sunlight,
- far from air conditioning systems,
- in a dust-free ambient.

Ambient temperature min.: +10 °C^(*)
 max: +30 °C
Relative humidity RH max: 80 %
 (*) The ambient temperature must not be inferior to the chosen set-point.

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



A= Lato muletto
 B= Lato transpallet

GB 3 TECHNICAL SPECIFICATIONS

3.4 TECHNICAL LAY-OUT DIAGRAMS

A= Fork lift truck side
 B= Transpallet side

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GB 4 HANDLING AND REMOVAL OF PACKAGING

4.1 PERSONNEL REQUISITES

Personnel in charge of handling the machine need no special requisites (just remember the type of packaging).

However, we suggest that this is done by someone who is accustomed to using machines for lifting and transport operations.

4.2 MACHINE CONDITIONS

The machine is normally supplied packed and on a pallet. If the machine is delivered by our staff, it may be without packaging. Other types of packaging may be supplied according to the destination and/or customer's needs.

4.3 EQUIPMENT NEEDED FOR HANDLING

The machine can be raised and moved by a lift truck, bridge crane, crane or other suitable means with an adequate capacity.

Check that the forks are level and the load is stable before moving the machine.

In order to handle the machine, the following rules should be observed:

- Move the machine slowly
- Do not tilt the machine, always keep it in an upright position
- Always make sure you can control the machine during handling

Warning:

- Do not turn the machine upside-down
- Do not drag the machine
- Do not shake the machine


How to store the machine

it should be kept in a dry environment with an ambient temperature ranging from 0° ÷ 40° C.

Do not place packaged machines one on top of the other and always keep them in an upright position as indicated by the special arrows marked on the packaging itself.

Make sure that the forks do not hit the machine frame. If it is hit accidentally, check immediately for any damage and if necessary, contact the manufacturer.

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GB 4 HANDLING AND REMOVAL OF PACKAGING

4.4 HOW TO RAISE AND MOVE WITH A BRIDGE CRANE

Accessories: textile fibre belts and PVC angle bars.
Weight: see technical specifications

Diagram of slings

- 1 Without packaging
- 2 With packaging

Do not use place the machine in a sling of metal cables or chains as these could damage it. Never put a sling around a machine without its pallet. In order to prevent the oil in the compressor from flowing into the refrigerating circuit, only transport, stock and handle the equipment in an upright position according to the instructions on the packaging. If the machine is placed on its side, leave it in an upright position for at least 24 hours before switching on.









4.5 REMOVAL OF PACKAGING

The machine is delivered ready-to-go. Before removing the machine from its packaging, check that it has not been damaged during transport. If so:

- Contact your area representative
 - Make a written report and if possible attach photographs showing the damage
 - Send a copy of the report to:
 - Shipping agent insurance company
 - Shipping agent
 - Manufacturer or area representative
- Remove all the packaging and destroy according to the laws in force.
 - Eliminate the parts with nails and eliminate plastic bags and plastic sheeting.

All the packaging materials can be recycled and can be disposed of according to the regulations in force. Take care to eliminate the packaging so that it does not become a danger and throw away the plastic bags, as they could cause children to suffocate.

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GB 4 HANDLING AND REMOVAL OF PACKAGING

4.5.1 How to remove the outer straps (fig. 4.5)

4.5.2 How to remove the cardboard casing (fig. 4.6)

4.5.3 Removal of plastic materials (fig. 4.7)

 When removing the plastic materials, make sure that you don't scratch the surface of the appliance.

4.5.4 How to remove the wooden side blocks (fig. 4.8)

- Remove the screws that block the ledges on the cabinet.
- Use a special lever to remove the two wooden blocks.

 Take care not to scratch the surface of the machine.

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GB 4 HANDLING AND REMOVAL OF PACKAGING

4.5.5 How to remove the pallet

- Just raise the machine a few centimetres so that the pallet can be pulled out from underneath.
- The machine can be placed directly on the floor.

 **Do not bang or scratch the machine**
In case of accidental bumps, please check if there is any damage and, should it be necessary, contact the manufacturer.

MANUAL HANDLING

All models are supported by four wheels; it is therefore possible to move manually. Before moving the equipment:

- Bear in mind the mass which has to be moved and calculate

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the number of people needed to control the mass during movement.



- check that the floor is smooth and completely flat.

4.6 HOW TO BLOCK THE MACHINE

The machine must not be fixed to the floor. The machine has been designed to operate on a horizontal floor.


In order to prevent the machine from moving accidentally, use a spanner to level and stabilize the machine as follows:

- Remove the lower side panels L (from fig. 4.10 to 4.12). For the right panel, first open the plexiglas panel protecting the chart recorder (fig. 4.11).

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

GB 4 HANDLING AND REMOVAL OF PACKAGING

- Adjust the small front feet P (fig. 4.13 - 4.14) using a spanner (fig. 4.15).
- Reassemble the lower side panels and close the chart recorder plexiglas panel.

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

GB 5 DESCRIPTION OF THE SYSTEM

5.1 GENERAL VIEW





- 1 Access door
- 2 Chart recorder (optional)
- 3 Control panel
- 4 Condenser protection grid
- 5 Handle to open door
- 6 Door closed
 - 6.1 With the Bioguard® system (PLUS model)
 - 6.2 With a key
- 7 LN₂/CO₂ connection
- 8 Activation/deactivation of the Bioguard® system and locking/unlocking of the control panel (PLUS models only)
- 10 Supply cable
- 11 Spacers

- 12 Compensating valve (this compensates the difference in external and internal pressure and helps the door to open).
- 18 Main switch

The structure is completely metallic made of galvanised sheet steel painted with epoxy resins. The door 1 has the same features as the main body, inside 9 it is insulated with thermoformed, high impact polystyrene and has a magnetic gasket. The spacers 11, located behind the machine, ensure there is enough room between the freezer and the wall to allow a sufficient circulation of air.

- A Door closed with the BIOGUARD® system
B Door closed with a KEY

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GB 5 DESCRIPTION OF THE SYSTEM

5.2 INTERNAL VIEW

The inner storage compartment 15 is made of AISI 304 stainless steel.




The appliance is equipped with inner doors or lids 14 used to increase the thermal insulation grade.

- 1 Access door
- 9 Door internal side
- 14 Inner doors
- 15 Inner storage compartments
- 16 Shelves racks
- 17 Magnets

Warning! If keep open the inner door too long, the water drops are formed on the gasket. When the drops freeze, the ice can block the inner door, so clean and dry the gasket and the inner door knocker before close.

The shelves inside the main body can be adjusted in height according to the customer's needs.

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
GB 5 DESCRIPTION OF THE SYSTEM

5.3 DESCRIPTION OF EXTERNAL CONNECTIONS

- 7 LN₂/CO₂ connection
- 8 Activation/deactivation of the Bioguard® system and locking/unlocking of the control panel (PLUS models only)
- 10 Supply cable
- 13 Cable guide for supply cable (without box)
- 18 Main switch
- 19 Line fuse
- 20 Remote control alarm electrical connection box
- 21 RS232 / 422 / 485 connector (optional)

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GB 5 DESCRIPTION OF THE SYSTEM

5.4 CONTROL SYSTEM

The control panel consists of a microprocessor-run electronic board, divided into two sections: one to control temperature regulation and the other to run the alarm system and temperature visualisation.

Both parts run the refrigerating units.

This board has 5 relay outputs to pilot the mechanisms.


It is also equipped with three analog inputs for the PT100 probe (to control temperature, low stage condensation and probe alarm). An input has been prepared to connect a buffer battery (12 V, 2 Ah) in order to guarantee the alarm system for over 48 hours in a blackout.

How to identify the keys

- 1 Not used.
- 2 Appliance ON/OFF key with relative warning light (green)
- 3 LED CONDENSER: if the led warning light flashes, the condenser is blocked or the ambient temperature is equal or higher than +32 °C.
- 5 Temperature setting key (control)

- 6 Temperature increase key (control)
- 7 Temperature decrease key (control)
- 8 LED COLD: cooling warning light (green)
- 9 LED POWER: warning light to show mains switched on
- 10 Set temperature display
- 11 Internal temperature/alarm set display
- 12 Flashing led: minimum temperature alarm
- 13 Flashing led: any alarm condition
- 14 Flashing led: maximum temperature alarm
- 15 Flashing led: battery OFF
- 16 Flashing LED warning light: the door is open (this operates only for machines with a key lock)
- 17 Flashing led: lack of line voltage
- 18 Temperature increase alarm key
- 19 Alarm set/alarm ON-OFF hidden key
- 20 Temperature decrease alarm key
- 21 ON/OFF alarm key with relative warning light (green)
- 22 Key to switch off acoustic alarm with relative warning light (green) / Alarm reset in the memory
- 23 If this key is pressed, it starts the alarm test
- 24 Acoustic alarm (buzzer)

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GB 5 DESCRIPTION OF THE SYSTEM

5.5 BIOGUARD® SYSTEM (only PLUS models)

This is an innovative system which enables the door to be opened by means of a *transponder* with a personal card (see. par. 8.6 - 8.7).


The system can also operate using a 12 Vdc, battery when no mains supply is available to give a 24 hours autonomy.

- 1 Led warning light to signal door position.
 - 2 Buzzer to signal door open
 - 3 Sensitive area for door control
 - 4 Personal card (TAG)
n° 1 MASTER card (red)
n° 3 USER cards (blue)
- A Front side
B Rear side

5.6 COOLING SYSTEM

The cooling system is completely sealed and uses two particularly silent compressors.

To obtain low evaporation temperatures (-90 °C) a cascade cooling system is used; it consists of two separate systems using two different refrigerants.

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GB 5 DESCRIPTION OF THE SYSTEM

5.8 VOLTAGE BOOSTER (optional)

The freezer can be equipped with a voltage control system named "VOLTAGE BOOSTER".

The system consists of:

- a microprocessor electronic board to measure the mains voltage;
- a transformer capable of withstanding inrush current;
- a contactor.

If the voltage drops, the system is capable of returning voltage to its pre-set values.

At times of peak voltage, the electronic circuit is cut off in order to protect the machine which will then stop. When the mains voltage returns within the pre-set limits, the system is automatically reset.

5.9 FUNCTIONING OF THE AUXILIARY LN₂/CO₂ COOLING SYSTEM (optional)

The auxiliary LN₂/CO₂ cooling system is used to cool the internal refrigerating compartment down by spraying liquid LN₂/CO₂ through a capillary pipe; the expansion of this gas/refrigerant allows the freezer to reach very low temperatures (down to -75°C). If the freezer is supplied with the above system, whenever the internal temperature accidentally exceeds the alarm limit, an electrically-driven valve opens the circuit allowing the passage of LN₂/CO₂ into the freezer.

Liquid LN₂/CO₂ (carbon dioxide) being a throwaway material is used only when no other solution can be found to keep the stored products safely. The auxiliary LN₂/CO₂ cooling system is fed by an automatically recharged buffer battery which permits over 48 hours of operations, even during a black-out.

WARNING! If the buffer system has been working for over 48 hours, 12 hours are needed to recharge it in order for it to be efficient once again.

5.10 INSULATION

Mod. 340 - 500 - 110 - 370 - 550

Average thickness: 140 mm. Obtained by injection of polyurethane foam (40 kg/m³).

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GB 6 INSTALLATION

6.1 POSITIONING OF THE MACHINE


- Check the ambient conditions of the installation site (see paragr. "3.3 ENVIRONMENTAL CONDITIONS").
- In order to bring the machine to the installation site, refer to chapter 4.
- During installation remember to leave enough space for work and maintenance operations.
- Check that the floor is perfectly flat.

WARNING:

The machine has not been designed to operate in areas classed as dangerous (explosive), provided for in the CEI (EIC) norm 64/4. The machine cannot operate in an explosive or corrosive ambient. The equipment being classified as IP 20 is not fit for operation outdoors or in places that are not protected against atmospheric agents.

Should the user need to place the equipment in places other than those described above, always contact the manufacturer before starting operation.

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GB 6 INSTALLATION

6.2 HOW TO CONNECT THE LN₂/CO₂ CYLINDERS (optional)

- Set up the cylinder connection as shown in the diagram (showing all the component parts for correct installation). Use liquid LN₂/CO₂ cylinders with a ball cock valve with plummets.

The average consumption of a freezer whose alarm limit is set at -70 °C is 2,5 kg/hour approximately (for CO₂).


The cooling system is guaranteed to run for over 48 hours even during a black-out, thanks to a buffer battery, which automatically recharges when the mains supply returns.

- 1 Freezer
- 2 Connection hose pipe (supplied)
- 3 Mechanical filter (not supplied) (*)
- 4 Manual check valve (not supplied) (*)
- 5 Gauge (not supplied): it just indicates whether there is LN₂/CO₂ in the cylinders (*) It does not indicate the quantity, because the pressure is kept constant as long as even a minimum quantity of liquid is present.
- 6 Header for cylinder connection (not supplied) (*)
- 7 Cylinder/header connection hose pipe (not supplied) (*)
- 8 Liquid LN₂/CO₂ cylinders bottles complete with plummet (not supplied) (*)


* For customer's use



WARNING! IMPORTANT NOTE!
If the freezer is equipped with the auxiliary CO₂ cooling system, the installation site must be well aerated in order to avoid an excessive accumulation of carbon dioxide, which should be dangerous if inhaled in big quantities.

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

GB 6 INSTALLATION

To attach the CO₂ connection tube 2 to the cylinder(s):

- insert gasket G into adapter nut D;
- screw the adapter nut D on to the cylinder;
- put pipe neck C into position and screw the pipe union B to the adapter nut D.

To attach the LN₂ connection tube 2 to the freezer:

- unscrew the valve cap as shown in the figure;
- screw pipe union B 1/4 SAE GAS onto the LN₂ connection.

LN₂ valve operating limits:

- min 0 bar;
- max 2.0 bar.

The connection tube can withstand a pressure of 100 bars.

6.1.1 HOW TO REPLACE THE CYLINDER

When the cylinder runs out of LN₂/CO₂ and has to be replaced, check that the pipe neck C of the connection tube 2 (fig.6.5) is not damaged (fig. 6.6).

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Fig. 6.6 = Cartella INTEGRA
Fig. 6.7 = Cartella DETERIORATA

If the pipe neck appears damaged e.g. it shows cracks (fig. 6.7), excessive wear of the copper layer etc., replace the connection tube 2.

The connection tube should, in any case, be replaced after the LN₂/CO₂ cylinder has been replaced three times.

The connection tube can withstand a pressure of 100 bars.

In order to replace the cylinder connection tube 2:

- unscrew the pipe union B;
- unscrew the adapter nut D.

To order the connection tube 2, ask for "LN₂/CO₂ CONNECTION PIPE"; specify also the model and registration number of your equipment (ref. Chapter 1).

Fig. 6.6 = UNDAMAGED Pipe neck
Fig. 6.7 = DAMAGED Pipe Neck

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GB 6 INSTALLATION

6.3 ELECTRICAL WIRING

The appliance must be supplied with electricity with the same tension as the one specified on the data tag (see chapter 1 "SUMMARY OF RATING PLATE DATA").

Class 1 insulation, complete with supply cable, without differential switch protection degree \leq 0.5 mA. "THE EQUIPMENT COULD BE A SOURCE OF ELECTRIC DISCHARGE IN THE EVENT OF BREAKDOWN OR INTERRUPTION OF THE PROTECTION CONDUCTOR. THEREFORE, TO IMPROVE SAFETY CONDITIONS, IT COULD BE ALSO INSTALLED (ACCORDING TO THE USER'S REQUIREMENTS) PERMANENTLY" (ref. 6.12.2.1 of CEI/EN 61010-1 adjustments).

This equipment can be considered electrically safe only when it has been correctly wired and an efficient earth system has been installed as foreseen by the safety norms in force.

This fundamental point of safety must be checked and if you are in any doubt, ask for the system to be accurately checked by professional qualified personnel.

In order to ensure that the appliance has been correctly and safely installed, a special, easily accessible socket must be placed near the appliance that is controlled by a multi-polar switch with a 3 mm or more contact opening, in compliance with the safety norms in force, equipped with fuses, with associated differential.



Certain fundamental rules must be followed when using this machine as when using any electrical machine, such as:

- do not touch the equipment with wet or damp hands or feet;
- do not place the machine on normally wet floors;
- do not use extensions in bathrooms or showers;
- do not pull on the supply cable in order to disconnect it from the mains supply.

6.3.1 How to connect to the electrical mains supply

NOTE: Machine sectioning takes place by disconnecting the plug from the socket.

The machine is supplied with a cable already connected to the panel terminal board; just insert the plug into the mains socket.

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GB 6 INSTALLATION

6.3.2 Replacement of the supply cable

During this operation the equipment has to be disconnected from the electrical supply. In order to disconnect the equipment from the electrical supply extract the plug from the socket.

The supply cable 10 complete with 16A plug in compliance with the adjustments in force (CEI-EN 61010-1 ref. 6.12.2.2) is connected directly to the power container by means of a cable guide.

WARNING! It is strongly recommended to unroll the supply cable for all its length and check that it is not crushed.





The supply cable must be replaced by qualified personnel.

6.3.3 Remote alarm connection

The attachment 20 for the remote alarm is a clean switch contact. The contact can tolerate a maximum electricity of 2 A at 230 V.

The contact exchanges after 30 seconds under a maximum alarm or minimum temperature condition, and after 10 seconds when the mains power is switched off.

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
GB 6 INSTALLATION

6.3.4 Serial interface (optional)

3 types of connections are available: RS485 (default), RS232 (on request) and RS422 (on request). Unless otherwise specified, the chamber is supplied with an RS485 interface. In this case the chamber has only to be connected to the relevant serial port on the PC. All available lines are to be found on the connector. The connecting cable is supplied according to the configuration requested by the customer as shown in figures 6.11-12-13.

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GB 7 USE FORESEEN BY THE MANUFACTURER

7.1 AIM OF THE MACHINES

These machines are en-bloc refrigerating cupboards, used to conserve diagnostic drugs, biological material, serums and vaccines.

They are furthermore appropriate for the preservation of industrial products in general, within the temperature zones indicated in chapter "3 technical characteristics".

For further information please refer to chapter 3 entitled "Technical specifications".

7.2 INTENDED USE OF THE MACHINES

The machines are intended to be used mainly by research laboratories and should only treat the types of materials mentioned in paragraph 7.1.

Before treating any different types of materials, you should consult the manufacturer directly to obtain his specific authorisation.

7.3 OPERATOR

No special technical knowledge is needed to use this machine.

7.4 USE LIMITS

The machines can be used within the ambient temperature and humidity range described in chapter 3 entitled "TECHNICAL SPECIFICATIONS". This series of machines must not contain:


- inflammable and/or explosive liquid or solid materials,
- liquid or solid inflammable materials.
- liquid or solid materials which could generate inflammable gases during evaporation or sublimation.
- liquid or solid materials containing chemically aggressive substances,
- liquid or solid materials could generate aggressive fumes during evaporation or sublimation.

7.5 USE OF PROTECTIVE CLOTHING

Suitable protective clothing must be worn, as stated in the norms in force, for protection against the temperatures of both the inner storage compartment and of the frozen products themselves.

7.6 DANGEROUS AREAS AND WASTE RISKS


When the equipment is switched on, personnel must take into consideration the risks from low temperatures. Take all the necessary precautions and wear suitable protective clothing.

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the Bioguard® system and locks the keyboard). During machine storage, Bioguard® is switched off in order to prevent the buffer battery from discharging completely.

8.3 SETTING THE OPERATING TEMPERATURE In order to set the operative temperature:

- Press key 5 to see the operational temperature (fig. 8.3).
 - Use keys 6 and 7 (increase and decrease) to modify it. The value is shown on display 10.
- The display switches off after a few seconds. When this operation has been completed, the cupboard will return to the set temperature value.

GB 8 START-UP

8.2 START UP

At the time of installation:

- insert the plug into the main socket;
- turn the main switch 18 to ON;
- press key 2 (ON-OFF) on the control panel to turn on the machine;
- set the operating temperature and with the alarm system switched off wait for the instrument to reach the set temperature (see par. 8.3);
- (Optional) turn the key 8 to the  symbol position (activates the Bioguard® system and locks the keyboard). Then memorise the card as described in par. 8.6 (Bioguard® model only). When the machine has reached the set value, turn the key 8 to the  symbol position (deactivates the Bioguard® system and unlocks the keyboard).
- Press key 21 together with key 19 (on the control panel - alarm section) to activate the alarm (Bioguard® model only).
- Set the alarm temperature setting (see. par. 8.4).
- (Optional) turn the key 8 to the symbol  position (activates

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GB 8 START-UP

8.4 SETTING THE ALARM TEMPERATURE SET POINT

NB: the maximum and minimum temperature alarms act as a safety system for material being treated: we advise you to set the minimum temperature alarm to trigger at -5 °C below the set point and the maximum temperature alarm at +10 °C above the set point. The alarm system is switched on or off by pressing key 19 and key 21 at the same time. We advise you to set the alarms when the machine has reached its working temperature, otherwise the CO₂/LN₂ valve (optional) will trigger.

- Press key 19 (SET ALARM) once in order to set the minimum alarm; twice in order to set the maximum alarm. The respective LED warning lights will flash and the threshold value will appear on the display. Use keys 18 and 20 to modify the value. When this operation has been completed, the cupboard will return to the set temperature value.

WARNING! When the maximum temperature alarm is triggered, the CO₂/LN₂ solenoid valve (optional) also triggers.


In order to switch off the buzzer during an alarm condition, press key 22 (BUZZER OFF); the buzzer will continue to beep every 15 minutes if the alarm conditions have not ceased. The buzzer starts to work under all alarm conditions, with the exception of the blocked condenser alarm.

NOTE: when an alarm condition is detected, the visual signal is sent first and the acoustic signal follows after 30 seconds; the relay for remote alarm transmission is closed simultaneously.

At the end of any alarm condition, the corresponding led light remains switched on: in order to reset, press key 22 (buzzer).

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


GB 8 START-UP

8.5 ACTIVATION/DEACTIVATION OF THE BIOGUARD® SYSTEM AND LOCKING/UNLOCKING THE CONTROL PANEL (Only for models equipped with this device)

- Turn the key 8 to the symbol  position:
- the Bioguard® system is activated, allowing card memorization;
- the control panel is locked (the displaying of the set point values and the switching off of the buzzer can be done by pressing key 22 on the control panel).

	CONTROL PANEL	BIOGUARD® SYSTEM
	Key lock: inactive It is possible to: - use the keys on the control panel	DEACTIVATED It is possible to: - open the door in emergency conditions.
	Key lock: active It is only possible to: - display the set point values - switch off the buzzer with key 22	ACTIVATED It is possible to: - memorize cards.

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GB 8 START-UP

8.6 CARD MEMORISATION (TAG) - (Bioguard® option only)

be shown in order to exit from this process. Certain important precautions should be taken at this stage to guarantee the unit will operate correctly.

 NB: the system is supplied with 1 MASTER card + 3 USER cards. Extra USER cards can be purchased as an accessory.

- A TAG front**
Blue background: USER card, used to open the door.
Red background: MASTER card, used to memorise the USER card, *it does not open the door.*
- B TAG back:** space to write the name of the USER.



- Place the MASTER card in front of the reader. After the beep the led light turns orange.
- Place up to a maximum of 10 USER cards in front of the reader. A longer beep indicates that each card has been read correctly.

The system automatically exits from the memory mode when the tenth card authorised to open the door is placed in front of the reader. If there are fewer than 10 users, the MASTER card will have to

In particular:
- when all the cards have been memorised, make sure that **ALL** the USER cards have been correctly read and memorised by the system by placing each one in front of the unit door.

If they have been memorised incorrectly, the MASTER card must be replaced in front of the reader in order to cancel all the codes and all the cards must be passed once again in front of the reader. The same procedure must be followed if authorisation has to be denied to one of the cards. The system can memorise up to a maximum of 10 USER cards and the pre-set codes can be reset after a blackout. This procedure has to be followed when the machine is first started or each time further USER cards have to be memorised.

 **Warning! Do not use the master card as user card. This use erase the memory system.**

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GB 8 START-UP



8.7 HOW TO OPEN/CLOSE THE DOOR

8.7.1 Description of how to open/close the door with a key

Open the door as indicated in figures 8.79 through 8.12.
The led light 16 on the control panel (fig. 8.1) indicates when the door is open.
In order to close the door, follow the procedure in reverse.

8.7.2 Description of how to open/close the door with the BIOGUARD® system

When the door is closed, the led warning light 1 is switched off (fig. 8.13).
When the card is passed in front of the *Bioguard*®, an acoustic sound is heard and the led light 1 turns green (fig. 8.14).
You now have 10 seconds in which to open the door (fig. 8.15).
When the door is open, the led light turns red; after 30 seconds you will hear an acoustic signal.
When the door is closed, the red led light 1 flashes and the door cannot be opened for 15 seconds.

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









GB 9 USE

9.1 HOW TO INSTALL THE DRAWERS (optional) - vertical models only

- Insert the guide-runs A into the racks C (front and rear, on the right-hand and left-hand side) as shown in figs. 9, 9.1. In order to insert drawer B, proceed as follows:
- insert the drawer wheels RT behind the guide-run wheels RA, as shown in figure 9.3. This system prevents the drawers from sliding out of the guide-runs when they are opened.
- Slide the drawers B along the guide-runs A and push them inside the chamber.

- A Guide-runs
- B Drawers
- C Racks

⚠ WARNING! If the machine has been loaded (i.e. there is material in the drawers) open one drawer at a time. Do not open all the drawers at the same time, otherwise the excess weight may unbalance the machine and cause it to fall over face down.

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GB 9 USE

9.2 HOW TO INSTALL THE SHELVES - vertical models only

- Follow the operations shown in figs. 9.6 and 9.7 in order to mount supports S at the required height.
- Rest the shelves P on the supports S and check:
 - their stability;
 - that they are resting firmly on all four supports S.

9.3 PRODUCT LOADING

Insert a few objects at a time inside the refrigerating compartment after temperature has become stable. Ambient temperature material (e.g. +25 °C) should be put inside the freezer gradually and in small quantities to ensure proper storage; allow the equipment to freeze the products inside the refrigerating compartment before introducing new material to be frozen. If this important rule is not observed, the temperature inside the freezer could suddenly increase and cause previously stored products to deteriorate. Whenever material is taken from or put inside the freezer, act very quickly since door opening determines a continuous change of air with subsequent temperature increase and ice formation on internal walls.

Do not store products imbued with liquid substances inside the freezer since they could produce vapours during the first cooling phase, leading to ice formation on internal walls.

WARNING! Do not store liquid products in hermetically-sealed containers as freezing could cause the containers to break and the liquid to leak.

9.4 TEMPERATURE CHART RECORDER (optional)

Freezers can be supplied with a temperature chart recorder 2 upon customer request.

The standard 7-day battery-operated recorder comes complete with a clock run by a 1.5V high-capacity accumulator lasting approximately 6 months. The temperature chart recorder is moreover supplied complete with no. 1 printing cartridge and no. 55 chart disks ensuring 1-year coverage. Contact Angelantoni's Assistance Service for additional printing cartridges and chart disks.

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GB 9 USE

9.4.1 Use of the recorder


The recorder has a battery powered clock.

Chart disk replacement.

- Open the plexiglas door G.
- Lift the pen 4.
- Move the lock 1 in the direction shown by the arrow and lift it.
- Remove the chart disk and install the new chart disk in the central pin. Move the lock 1 to the previous position and lower the pen 4.

Battery substitution (every 4 months) (fig. 9.16).

- Follow the steps described in the figs. 9.13 - 9.15, and move the diagram disk without removing it.
- Replace the battery when necessary with a new 1.5 V (Mod. AA LR6) battery.

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GB 9 USE

WRITING
Writing cartridges.
 The instruments are supplied with pen points already inserted into their respective rods. In order to start up, proceed as follows:


- remove the protective, plastic caps from the cartridge writing tips.

Writing cartridges replacement
 If the writing should appear faulty, replace the pen points. Proceed as follows:

- Hold the rod firmly and remove the pen point that has to be replaced by pulling it downwards.
- Insert a new pen point of the same colour and the same length.

To order the nib, please specify the code as indicated in the following table.
 The code number to order the recorder disk is printed on the centre of the recorder disk itself.

RECORDER TYPE	CARTRIDGES	
	FIRST PEN	SECOND PEN
SINGLE PEN	RED MEDIUM POINTS	
TWO PENS	RED MEDIUM POINTS	BLUE LONG POINTS

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GB 9 USE

9.5 DATA DOWNLOADING VIA SERIAL PORT (optional)

- Connect your computer to the serial port 21 (see paragr. 6.3.5).
- Switch on your computer and start the programme.
- In order to set and start up the software, please refer to the attached instruction handbook according to the configuration of your machine.

A : Cryolog™ software instruction handbook

9.6 Additional probe (optional)
 It's possible to buy an additional temperature probe. This probe could be connected to an external device and allow to have the temperature measurements. The probe could be supplied, if requested, with a standard connector.

9.7 Port hole (optional)
 The freezer could be supplied with a port hole (fig. 9.19) that allow the entry of external temperature measurement devices.


9.8 WATER CONDENSER CONNECTION (optional)
WARNING! This description does not concern machines with an air condenser. The connections for water condenser are located in the rear part of the equipment.

- Connect the pipes to the connections referring to the label. This chamber is built only for the choosed type of water when purchasing. The minimum water flow rate must be 0,4m³/h.

WATER CHARACTERISTICS

Pressure bar (1bar: Pa x 10 ⁵)		Temperature °C	
min	max	min	max
2	6	10	20

Regularly check the hydraulic connections and pay special attention to the hose pipes.
 The connection pipes between the water plug and the chamber must be hose pipes in order to allow partial movement of the chamber itself during maintenance. Installation is advisable for regular cheks of: thermometer, pressure gauge with pressure reducer and manual ON/OFF valves.

A	110	370	550	340	500
380 / 7 / 50	9				
380 / 5 / 75	9				
600 / 10 / 50		21	33		
600 / 7 / 75		21	33		
171 / 12 / 50				18	30
161 / 8 / 75				18	30
401 / 28 / 50				9	15
404 / 20 / 75				9	15

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GB 9 USE

9.9 STORAGE CONTAINERS

The racks were designed to conveniently organize the storage of products inside the chamber.

A sample unit: 161 / 8 / 75

- 161 Container height;
- 8 Number of boxes container in a single container
- 75 Height of a single box

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
GB 10 SAFETY DEVICES- CHECK AND SET-UP

SAFETY DEVICE	CHECK	SET-UP
Key to lock handle (B)	This key ensures that the machine cannot be opened either accidentally or on purpose by unauthorized personnel. The chamber may be at very low temperatures.	<ul style="list-style-type: none"> • Close the door, remove the key and check that the door will not open.
Electromagnetic lock (Bioguard® system) (A) - Optional	This closing ensures that the machine cannot be opened either accidentally or on purpose by unauthorized personnel. The chamber may be at very low temperatures.	<ul style="list-style-type: none"> • Close the door and check that the door will not open.

 The safety devices systems concerning the cooling system cannot be checked or set by the user and are, therefore, described in chapters 11 and 12. There are no more safety devices in the present machines.

- A Door closed with the BIOGUARD® system**
- B Door closed with a KEY**

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GB 11 MAINTENANCE

Maintenance must be carried out by qualified personnel.

Some maintenance operations may be carried out when the machine is working and therefore all the necessary safety precautions should be taken. With the exception of specific cases maintenance operations must be carried out in the following conditions:

- machine switched off
- chamber temperature the same as ambient temperature
- machine disconnected from the mains supply

11.1 HOW TO DEFROST AND CLEAN THE MAIN BODY AND DOORS

This operation should be performed every 3 months at least. This interval may however vary according to the number of times the door is opened in that period of time.

- Turn the general switch onto OFF (refer to 18).
- Disconnect the plug from the mains socket
- Open the external door and the inner doors of the freezer. After 2 hours, the freezer should be completely defrosted. Collect the water at the bottom of the freezer.

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


WARNING! Never spray water directly on any part of the equipment since the electric insulation could be seriously damaged and cause malfunctioning. Never try to remove the ice with your hands. Danger of scratches and burns from very low temperatures.

- Remove the internal shelves.
- Clean the main body, the door and the inner doors with a damp cloth dipped in a neutral detergent and then rinse with a damp cloth dipped in water (dirt deposits on door and inner door gaskets may compromise their functionality).
- Clean the Bioguard® magnet, if present, to avoid malfunctioning.
- Introduce the shelves again inside the freezer.
- Start the freezer again and store products only after the set temperature has been reached.

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GB 11 MAINTENANCE

11.2 CLEANING OF THE AIR CONDENSER

Clean the condenser every 2 months with soft brushes or blowing compressed air on it.

In order to carry out this operation proceed as follows:

- remove the screw V;
- remove grid 4 by pulling it toward you, using a screwdriver or other non-abrasive tool as shown in the photo 11.2.1-11.3. Be careful not to drop the grid.
- remove the inner filter (if present) and clean it;
- clean the condenser (lying just behind the grid);
- replace the grid.

11.3 CLEANING OF THE WATER CONDENSER

The water condenser cleaning (only for equipment that use this device) must be done only by qualified personnel.

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11.4 COOLING SYSTEM MAINTENANCE

Each machine is equipped with a refrigerating system using a hermetically sealed compressor, and therefore, (if there are any accidental leaks) no refrigerant or non-freezing oil has to be added. We advise you to have it checked once a year by qualified personnel.

11.5 MAINTENANCE OF THE ELECTRICAL-ELECTRONIC PART

The electronic control is a very simple, reliable, integrated system in a single container, so no maintenance is required.

We advise you to have it checked once a year by qualified personnel in order to see that the alarm system is operating correctly, that the visualised temperature remains linear, and to replace the compressor control relay (every 5 years).

Our company provides its own economic, subscriber "SPECIALIZED TECHNICAL ASSISTANCE SERVICE".

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
GB 12 TROUBLESHOOTING

The maintenance operations authorizes are those indicated by the producer are indicated in the "Regular maintenance" chapter. The information given in the column "SOLUTION" in the tables below does not authorize any operations if these compromise safety; the information is given to help any specialized technicians find the fault.


PROBLEM	PROBABLE REASON	SOLUTION
The machine does not start	<ul style="list-style-type: none"> • There is no supply from the mains. • The main switch is in "OFF" position • The plug is not connected to the main • Some other reason not mentioned above 	<ul style="list-style-type: none"> • Check and proceed accordingly • Turn ON the main switch • Connect the plug to the main • Call the technical assistance service.
The equipment does not reach the set temperature	<ul style="list-style-type: none"> • Setpoint programming error. • There is too much material in the test chamber. • The test compartment contain too wet material. • Malfunctioning of the regulation system. • Malfunctioning of the cooling system. • Condenser clogged 	<ul style="list-style-type: none"> • Check the set temperature value following the procedure as in (ref. par.8.1). • Reduce the quantity. • Allow temperature values to stabilise. • Call the technical assistance service. • Call the technical assistance service. • Clean the condenser; • Call the technical assistance service.
The equipment is too noisy	<ul style="list-style-type: none"> • The equipment is not evenly laid on the floor. • Contact with foreign bodies. • Damaged mechanical parts. 	<ul style="list-style-type: none"> • Remove the cause. • Check that foreign bodies such as boxes or other structures do not come into contact with the equipment so as to avoid vibrations during functioning. • Check that moving mechanical parts such as fans are not sources of noise.
The display shows alarm codes	<ul style="list-style-type: none"> • Refer to paragraph (ref. 8.1) of these operating instructions. 	<ul style="list-style-type: none"> • Call the technical assistance service and specify the type of code displayed as well as the possible cause.

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GB 12 TROUBLESHOOTING

PROBLEM	PROBABLE REASON	SOLUTION
The door does not open.	<ul style="list-style-type: none"> • Card has not been memorised or reading is incorrect. • Disabled code. • The card does not work. 	<ul style="list-style-type: none"> • Repeat the procedure to memorise the card (see paragr. 8.2). If the door remains shut, call your technical assistance service. • Repeat the procedure to memorise the card (see paragr. 8.2). If the door remains shut, call your technical assistance service. • Turn the key 8 to the  position. If the door remains shut, call your technical assistance service.
Bioguard® not receiving power.	<ul style="list-style-type: none"> • Power supply out of order. • Defective buffer battery. 	<ul style="list-style-type: none"> • Call the technical assistance service. • Call the technical assistance service.

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GB 13 REMOVAL FROM INSTALLATION SITE

The machine can be disassembled in order to move it from one site to another or in order to be scrapped. Transport must be carried out in compliance with the norms described in this handbook in the chapters "HANDLING" and "INSTALLATION".

13.1 DISASSEMBLY

- In order to disconnect the machine from the mains supply, follow the procedure described for the electrical wiring in reverse.
- Clean the surfaces accurately, oil the parts that could rust and cover the machine so as to protect it from dust, dirt and humidity.

13.2 SCRAPPING

- In order to disconnect the machine from the mains supply, follow the procedure described for the electrical wiring in reverse.
- Discharge all the refrigerant from the refrigerating circuit.
- Remove all the oil from the motors.


WARNING! This product should not be left in the surrounding environment:

- refrigerating substances;
- the compressor oil;
- polyurethane.

It would be advisable to contact specialised companies for the disposal and the recycling of the above materials.

- Divide the remaining parts according to their type and destroy in compliance with the laws in force.
- To remove and transport the machine please refer to the procedure described in the chapters "HANDLING" and "INSTALLATION".

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

GB 14 DIAGRAMS

14.2 REFRIGERATOR DRAWING (-85°C)

- | | | | |
|----|------------------------------------|----|------------------------|
| C | Door frame | 7 | Capillary high stage |
| E | Evaporator | 10 | Oil separator |
| H | High stage compressor | 12 | Low stage drier filter |
| L | Low stage compressor | 13 | Low stage capillary |
| M | Condenser ventilation motor | 14 | Expansion tank |
| P | Pre-condenser low stage | 15 | Capillary |
| P1 | Maximum pressure switch high stage | 16 | Intermediate exchanger |
| P2 | Maximum pressure switch low stage | | |
| 2 | Pressure control valves | | |
| 3 | Condenser | | |
| 5 | Propeller blade | | |
| 6 | Drier filter High stage | | |

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
GB 14 DIAGRAMS

14.3 ELECTRIC DIAGRAM (-40 °C) pag.1 of 2

- | | | | |
|----|-----------------------------|----|-----------------------------------|
| A | Electric board | S3 | Open door micro switch |
| A2 | Antinoise filter | T1 | Winter temperature switch |
| K1 | Compressor relay high stage | T2 | Obstructed condenser temp. switch |
| F1 | Line fuses | V1 | Line solenoid valve |
| G1 | Buffer battery | Z1 | Regulation probe |
| M1 | High stage compressor | Z3 | Alarm probe |
| M3 | Condenser fan | | |
| P1 | Max pressure switch | | |
| S1 | Main switch | | |
| S2 | Setting keys block switch | | |

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

GB 14 DIAGRAMS

14.4 ELECTRIC DIAGRAM (-85 °C) pag. 1 of 2

- A Control board
 - A2 Antijamming filter
 - F1 Line fuses
 - G1 Buffer battery
 - K1 High stage compressor relay
 - K2 Low stage compressor relay
 - M1 High stage compressor
 - M2 Low stage compressor
 - M3 Condenser fan
 - P1 High stage pressure switch
 - P2 Low stage pressure switch max
 - S1 Main switch
 - S2 Setting keys block switch⁽¹⁾
 - S3 Open door micro switch⁽¹⁾
 - T2 Obstructed condenser temp. switch
 - Z1 Regulation probe
 - Z2 Intermediate exchanger probe
 - Z3 Alarm probe
- ⁽¹⁾ Optional.

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
GB 14 DIAGRAMS

14.5 ELECTRIC DIAGRAM (equipment at -85 °C with BOOSTER) pag. 1 of 2

- A Control board E2003
- A2 Antijamming filter
- K1 High stage compressor relay
- K2 Low stage compressor relay
- F1 Line fuses
- G1 Buffer battery
- M1 High stage compressor relay
- M2 Low stage compressor motor
- M3 Condenser fan motor
- P1 High stage compressor pressure switch
- P2 Low stage compressor pressure switch max
- S1 Main switch
- S2 Setting keys block switch
- S3 Open door micro switch
- T1 Winter temperature switch
- T2 Obstructed condenser temp. switch
- Z1 Regulation probe
- Z2 Intermediate exchanger probe
- Z3 Alarm probe
- b Blue cable
- n Black cable
- w White cable
- r Red cable

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
GB 14 DIAGRAMS

14.5 ELECTRIC DIAGRAM (equipment at -85 °C with BOOSTER) pag. 2 of 2

- B CO2 feeder
- B1 Bioguard®
- C Bioguard® feeder
- G2 CO2 buffer battery
- G3 Bioguard® system buffer battery
- S2 ON/OFF feeder- ON/OFF Bioguard®
- V1 CO2 valve

SE3153_2_85_040505

DET NORSKE VERITAS

QUALITY MANAGEMENT SYSTEM CERTIFICATE

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Management Representative
 Vittore Marraugon

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DET NORSKE VERITAS

ENVIRONMENTAL MANAGEMENT SYSTEM CERTIFICATE

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Certificazione rilasciata in conformità al Regolamento Tecnico SINCERT RT - 09
 Certification has been granted in conformity with the SINCERT Technical Regulation RT - 09

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La validità del presente certificato è subordinata a sorveglianza periodica (ogni 6, 9 o 12 mesi) e al riesame completo del sistema con periodicità triennale (The validity of this certificate is subject to periodic audits (every 6, 9 or 12 months) and the complete re-assessment of the system every three years. Le aziende in possesso di un certificato valida sono promettee della ISO 9001 e della ISO 14001, tutti i diritti sono riservati. All rights reserved. Tutti i diritti sono riservati.)
ANGELANTONI INDUSTRIE S.p.A.
Det Norske Veritas Italia S.r.l.